

Amazing **All in one**
Facial treatment



Powerful performance / Mobility Compact Size / Safe and effective












Versatile Applications

- ✓ Deep Cleanses
- ✓ Face Renewal, Collagen Remodeling
- ✓ Anti Ageing
- ✓ Elimination of wrinkle
- ✓ Face rejuvenation



Integrated Technology

 AQUA	<p>Deep Cleanses & Exfoliates & Digital Peeling Deeply Infuses & Hydrates & Renews the skin Acne, Anti-Bacteria, Anti-Inflammation, Soothing effect</p>	 TOUCH	<p>7inch Color Touch LCD Convenient GUI</p>
 RF	<p>Collagen remodeling & Improve skin elasticity Safe and painless, with immediate effect</p>	 FOOT SWITCH	<p>Foot Switch (Convenient use)</p>
 ION	<p>Stimulates skin tissue, hair roots Gives skin elasticity & anti-aging</p>	 AUTO	<p>AquaPeel Automatic cleansing and power off function</p>
 ULTRASOUND	<p>Melts fat & cellulite cell, gives Smooth skin texture Helps to produce Collagen & Elastin</p>	 DIGITAL	<p>Flow and air control through digital mechanism Precise and convenient control</p>
 LED	<p>Blue LED 430nm</p>		

perfect **FACIAL CARE** **TOTAL solution**



Treatment Applications



RF



AquaPeel



Ion Lifting



Ultrasound



Mask Lifting



Mask Ultrasound





- Convenient drain box separation and cleaning
- Structurally durable platform



- Additional one touch solution
- Solution conservation
- Adjustment of patented technique



- Digital control solution
- Precise control solution (9steps)
- Reduces solution loss



- Automatic cleaning function
- Automatic power off function



- Detachable one touch connector
- Automatic program upgrade through USB



Specification

Tornado tip	Large, Small
Vacuum Range	Max 680 mmHg
Radio Frequency	1Mhz ,Multi Polar(6Pol)
Ion Lifting	500Hz (Digital Ion Lifting)
Ultrasound	1Mhz / 2W/cm ²
User Interface	7inch Color Touch LCD
Input Voltage	100~120V, 50/60HZ / 220~240V, 50/60Hz
Noise Level	45dB
Weight	12Kg