



6600 Vet Advanced

Floor-mounted Digital X-ray System for Vet

With 32kw high frequency generator, 6600 vet advanced is the economical DR system to meet clinical needs for vet applications.



WEBSITE

Specification

X-ray Generator

Output power	32kW
Tube voltage range	40~150kV
Current output range	10~400mA
Time of exposure	1ms~10000ms
Power requirement	220VAC

X-ray Tube

Focal spot	0.6mm/1.2mm
Focal power	20kW/50kW
Anode heat capacity	300kHU
Rotary speed	2800RPM

Flat Panel Detector

FPD type	a-Si/CsI
Image area	43×43cm
Acquisition matrix	2912×2912
Pixels size	139μm
Resolution	≥3.4lp/mm
A / D conversion	16bit

Integrated Table


Size(L x W x H)	1500×700×780mm
Horizontal moving range	400mm
Longitudinal moving range	100mm
SID	1.0m

Collimator


Automatic shutdown time	Provide
Light field	≥43cm X 43cm
Light resource	LED

Professional Workstation

User-friendly Interface



Full Measurement Tools



Display and operate in shared screen



Shenzhen Lanmage Medical Technology Co., Ltd.

No.103, Baguang Service Center, No.2 Baisha Bay Road, Baguang Community, Kuichong Subdistrict, Dapeng New District, Shenzhen, 518119, Guangdong, P.R. China

Website: en.lanmage.com
E-mail: oversea@lanmage.com



LINKEDIN



FACEBOOK