


**BOUNDLESS VISION
ENDLESS CARE**



ASR-4000E Digital Mammography System

Shenzhen Anke High-tech Co., Ltd.

Address: Block B, LingYa Industrial Zone, Tangtou No.1 Road,
Bao'an District, Shenzhen, 518108, P.R. China

Tel: +86-755-21622518 26688889

Fax: +86-755-26695307 26685908

Email: anke@anke.com

Website: www.anke.com

Follow us on:

-  ANKE(@ankeint)
-  ANKE(@ankemri)
-  ANKE(@ankemarcom)
-  Shenzhen Anke High-tech Co., Ltd





Super Configuration

High-precision Multi-function Rack

High accuracy movement control; x-ray tube rotates around detector isocenter; it provides hardware support for breast tomosynthesis and breast biopsy



The Best Choice for High Quality Mammography

The third generation full size amorphous Silicon (a-Si) flat panel detector with high x-ray conversion efficiency, better image resolution and contrast




Fully Automatic Exposure Technology

Precise and low dosage ensures the clarity of breasts



Excellent Image Quality




Innovative Human-care Design




Safety

Threshold pressure control technology ensures safety



Comfort

Multi-frequency self-induced compression technology ensures the best image quality and reduces patients' pain maximally



Easy Operation

High intelligent positioning by one button provides easy operation and improves work efficiency

